

HPLC Column

Affinity HPLC Column Selection and Handling Precautions

The choice of HPLC columns

The choice of column tubes

The materials of HPLC column tubes include SS316 (stainless steel), PTFE, and PMMA etc., determined by characteristics of mobile phase, pressure degree of column and sample. SS316 is used when mobile phase is organic solvents with pressure between 5 to 30 Mpa. When the mobile phase is 100% water or buffer solution with pressure less than 4 Mpa, PMMA or PTFE is chosen for less impact on activity of biological sample.

Inner diameter

A. 1-2mm ID column, is specific used for micro LC(MLC), such as LC -Mass spectrometry. But for routine analysis, it's not easy to use. Although the solvent consumption is small, the requirement is too rigid. It need the instrument only have a very small dead volume. Besides, this kind of columns is short-life.

B. 4-6mm (3.9, 4.0, 4.6, 5.0, 6.0mm) ID columns are analytical scale and suitable for routine analysis. 4.6mm ID columns are the most usual type. Best flow rate is 1ml/min which general instruments can match with. They have high column efficiency, stable performance, and longer life time.

C. 7.8-10.0mm ID columns are semi-preparative column. Chromatographic conditions can be transplanted from analytical column. They can be equipped on normal LC instruments to collect small amount of high purity components to quality and research.

D. 20-100mm ID columns are preparative column, which can prepare a large number of pure components with commercial value. At present, although the price is higher, it is must equipped for the pharmaceutical industry.

The length of the column tube

Length of HPLC columns is between 50 and 500mm. For general analysis 150-250mm is most commonly used. Columns longer than 250mm though have high column efficiency, have much higher pressure. So it is not economic just for better efficiency to increase the column length.

The choice of packing

Particle size

particle size of commonly used packing is 3 - 10um. Small particle size can achieve high column efficiency, but column pressure is also high. Column pressure is an important factor that cannot be ignored. High column pressure may lead to packing collapse and reduce columns' life time. Especially when mobile phase is methanol with larger water content, hydrogen reaction formed between water and methanol makes the viscosity to increase. If high column efficiency is pursued, water-acetonitrile system is recommended.

For preparative column, main pursuit is preparation volume, and separation is secondary. Generally packing with larger size than 10 um is chosen with low cost and low column pressure.

For UHPLC which has higher column efficiency, better separation and shorter separation time, particle size is so small that the pressure is much higher than HPLC. We offer two specifications: 1.8µm and 2.2µm. UHPLC columns can withstand pressures up to 10000psi, and have good reproducibility.

The specifications of packings

Molecular weight smaller than 2000

Aqueous samples	Non-ionic	Reverse phase chromatography	C30,C18,C18-WP, C8
	Ionic	Reverse phase chromatography	C30,C18,C18-WP, C8
		Ion exchange chromatography	SAX,SCX; Transgenomic ICSep AN, ICSep CN
	Amino acids	Reverse phase chromatography	C18
		Aqueous samples	Transgenomic AMINOSep amino acid column
	Organic acids	Ion exclusion chromatography	Transgenomic ICSep organic acid column;Shodex SUGAR SH1821, KC-811
	Monosaccharides, disaccharides, oligosaccharides, Peptide	Reverse phase chromatography	NH2
		Ion exclusion chromatography	Shodex SUGAR SH1821
Inverting chromatography		Transgenomic CARBOSep resin type sugar column	
Oil-soluble	Non-polar	Reverse phase chromatography	C18
		Reversed-phase chromatography	C18
	Polar	Normal phase chromatography	NH2, CN, SIL
Chiral sample		Normal phase chromatography	NH2, CN, SIL
		Chiral chromatography	Regis Whelk-O, RegisPack, RegisCell

CHROMATOGRAPHY

Molecular weight smaller than 2000

Aqueous samples	Non-ionic	Reversed-phase chromatography	C18-BIO
	Proteins, polypeptides	Gel filtration chromatography (GFC)	Shodex KW-800, SB-800 HQ; Gel X series
		Ion-exchange chromatography	SAX, SCX; Sep series
		Reversed-phase chromatography	C18-BIO
		Affinity chromatography	Shodex AFpak,
	Nucleic acid	Ion-exchange chromatography	SAX, SCX; Sep series
	Polysaccharide	Gel filtration chromatography (GFC)	Shodex KW-800, SB-800 HQ; Gel X series
Ion-exchange chromatography		SAX, SCX; Sep series	
Oil-soluble Samples		Gel filtration chromatography (GFC)	Shodex SB-800 HQ; Gel-S Series
		Reversed-phase chromatography	C18-BIO

Precautions for use of columns

Column Equilibration

When preparing to introduce your desired mobile phase into a new column, be aware of the miscibility of the solvents being introduced to the column and the solvent inside the column. If they are not, it is necessary to pump one or more miscible intermediate solvents through the column to avoid high pressure. Equilibrate the column with a minimum of 10 column volumes of mobile phase to be used.

Reversed-phase columns equilibration method

Reversed-phase columns equilibrate in as little as 20 column volumes of mobile phase. If the new eluent being introduced contains buffer salts, it is recommended that the column is flushed with a highly aqueous eluent (such as 90:10 Water: MeCN) before introducing buffer, to avoid precipitation of salts on the column. For extra precaution, introduce new buffered eluents WITHOUT the buffer component for 5-10 column volumes, and then switch to the fully buffered eluent composition. Precipitation of buffer salts on the columns is essentially irreversible and destroys the column. When switching between solvents with vastly different polarities, it may be necessary to first purge the column with a mutually miscible solvent such as Isopropyl Alcohol or Dioxane at a reduced flow rate (approximately 50% of normal). Flushing with a minimum of 5 column volumes is recommended (e.g. 10mL for a 150 x 4.6mm I.D. column).

Normal phase column equilibration method

Normal-phase columns require longer equilibration times (at least 50 column volumes). To ensure good reproducibility and faster equilibration of normal-phase columns, a small, constant percentage of water can be added to the mobile phase.

Column maintenance

Eluent pH: At pH above 8, silica gels begin to dissolve; at acidic pH below 2.0 certain bonded phases (particularly CN) become hydrolyzed and gradual loss of bonded phase can occur. While many customers use the columns outside both sides of the pH spectrum with excellent results and good column lifetime, the best lifetimes are usually obtained at intermediate pH conditions.

Pressure: To maximize column life operate at pressures up to 20 MPa (~ 3000 psi) for standard HPLC phases (UHPLC columns can be used at higher pressures, as indicated on the test chromatogram).

Sample Dissolution: Samples should be dissolved in the eluent or solvent weaker than the eluent, which helps avoid sample precipitation at the column head and inconsistent retention values. Filter sample with 0.45µm membrane to remove particulate matter before injection.

Solvents: Use HPLC or spectroscopy grade solvents that have been filtered through a 0.45µm filter. Filter all buffer solutions before use. Avoid introduction of particulates onto the column at all costs.

Guard Columns: Use a guard column of matching chemistry and particle size between the injector and main column. Guard columns need to be replaced at regular intervals as determined by sample contamination. When system backpressure limit, it is usually an indication that the guard column should be replaced. A sudden appearance of split peaks is also indicative of a need to replace the guard column.

Clean of Columns

Clean of reverse phase silica bonded phase columns

20 column volumes should be used for each wash stage:

95:5 water: ACN(Removal of buffer) → 100% ACN → 50:50 water: ACN

Clean of normal phase silica bonded phase column

20 column volumes should be used for each wash stage:

THF → Chloroform → Methylene Chloride → Hexane

Common Troubleshooting

Problem	Possible cause	Solution
No peaks or very small peaks	Detector off	Check detector
	Broken connections to recorder	Check connections
	No sample/Wrong sample	Check sample. Be sure it is not deteriorated. Check for bubbles in the vials
	Wrong settings on recorder or detector	Check attenuation. Check gain
No Flow	Pump off	Start Pump
	Flow interrupted	"Check reservoirs. Check position of the inlet tubing. Check loop for obstruction or air. Check degassing of mobile phase. Check compatibility of the mobile phase components."
	Leak	Check fittings. Check pump for leaks and precipitates. Check pump seals.
	Air trapped in the system	Disconnect column and prime pump. Flush system with 100% methanol or isopropanol. Contact servicing if necessary.
Column end leaks	Loose fitting	Tighten or replace fitting
	White powder at loose fitting	Cut tubing and replace ferrule; disassemble fitting, rinse and reassemble.
Leak at detector	Detector-seal failure	Replace detector seal or gaskets.
Leak at injection valve	Worn or scratched valve rotor	Replace valve rotor
Leak at pump	Pump seal failure	Replace pump seal; check piston for scratches and, if necessary, replace
Changing Retention Times	Buffer retention times	Use buffer with concentration greater than 20 mM.
	Contamination buildup	Flush column occasionally with strong solvent
	Equilibration time insufficient for gradient run or changes in isocratic mobile phase	Pass at least 10 column volumes through the column for gradient regeneration or after solvent changes
	First few injections - active sites	Condition column by injecting concentrated sample
	Inconsistent on-line mobile-phase mixing	Ensure gradient system is delivering a constant composition; compare with manually prepared mobile phase; partially premix mobile phase
	Selective evaporation of mobile-phase component	Cover solvent reservoirs; use less-vigorous helium purging; prepare fresh mobile phase
	Varying column temperature	Thermostat or insulate column; ensure laboratory temperature is constant.
	Active sites on column packing	Use mobile-phase modifier, competing base (basic compounds), or increase buffer strength; use higher coverage column packing.
	Column overloaded with sample	Decrease sample amount or use larger-diameter column.
	Increasing flow rate	Check and reset pump flow rate.
	Loss of bonded stationary phase or base silica	Use mobile-phase pH between pH 2 and pH 8
Increasing Retention Times	Decreasing flow rate	Check and reset pump flow rate; check for pump cavitation; check for leaking pump seals and other leaks in system
	Changing mobile-phase composition	Cover solvent reservoirs; ensure that gradient system is delivering correct composition.
	Loss of bonded stationary phase	Use mobile-phase pH between pH 2 and pH 8
Slow column equilibration time	Reversed phase ion pairing - long chain ion pairing reagents require longer equilibration time	Use ion-pairing reagent with shorter alkyl chain length
Void Time noise	Air bubbles in mobile phase	Degas or use back pressure restrictor on detector
	Positive-negative - difference in refractive index of injection solvent and mobile phase	Normal with many samples; use mobile phase as sample solvent
Drifting baseline	Negative direction (gradient elution) - absorbance of mobile-phase A	Use non-UV absorbing mobile phase solvents; use HPLC grade mobile phase solvents; add UV absorbing compound to mobile phase B.
	Positive direction (gradient elution) - absorbance of mobile phase B	Use higher UV absorbance detector wavelength; use non-UV absorbing mobile phase solvents; use HPLC grade mobile phase solvents; add UV absorbing compound to mobile phase A.
	Positive direction - contamination buildup and elution	Flush column with strong solvent; clean up sample; use HPLC grade solvents
	Wavy or undulating - temperature changes in room	Monitor and control changes in room temperature; insulate column or use column oven; cover refractive index detector and keep it out of air currents.
Baseline noise	Continuous - detector lamp problem or dirty cell	Replace UV lamp(each should last 2000 h; clean and flush flow cell.
	Gradient or isocratic proportioning- lack of solvent mixing	Use proper mixing device; check proportioning precision by spiking one solvent with UV absorbing compound and monitor UV absorbance detector output.
	Gradient or isocratic proportioning - malfunctioning proportioning valves	Clean or replace proportioning precision valves; partially remix solvents.
	Occasional sharp spikes - external electrical interference	Use voltage stabilizer for LC system; use independent electrical circuit.
	Periodic - pump pulses	Service or replace pulse damper; purge air from pump; clean or replace check valves.
	Random - contamination buildup	Flush column with strong solvent; clean up sample; use HPLC grade solvent
	Spikes - bubble in detector	Degas mobile phase; use back pressure restrictor at detector outlet.
	Spikes - column temperature higher than boiling point of solvent	Use lower column temperature.

Common Troubleshooting

Problem	Possible cause	Solution
Decreasing Pressure	Insufficient flow from pump	Loosen cap on mobile phase reservoir
	Leak in hydraulic lines from pump to column	Tighten or replace fittings; tighten rotor in injection valve
	Leaking pump check valve or seals	Replace or clean check valves; replace pump seals.
	Pump cavitation	Degas solvent; check for obstruction in line from solvent reservoir to pump; replace inlet-line frit
Fluctuating Pressure	Bubble in pump	Degas solvent; purge solvent with helium
	Leaking pump check valve or seals	Replace or clean check valves; replace pump seals
High Back Pressure	Column blocked with irreversibly adsorbed sample	Improve sample cleanup; use guard column; reverse-flush column with strong solvent to dissolve blockage
	"Column particle size too small (for example 3 micrometers)"	Use larger particle size (for example 5 micrometer)
	Microbial growth on column	"Use at least 10% organic modifier in mobile phase; use fresh buffer daily; add 0.02% sodium azide to aqueous mobile phase; store column in at least 25% organic solvent without buffer"
	Mobile phase viscosity too high	Use lower viscosity solvents or higher temperature
	Plugged frit in in-line filter or guard column	Replace frit or guard column
	Plugged inlet frit	Replace endfritting or frit assembly
	Polymetric columns - solvent change causes swelling of packing	Use correct solvent with column; change to proper solvent composition; consult manufacturer's solvent-compatibility chart; use a column with a higher percentage of cross-linking
	Salt precipitation (especially in reversed-phase chromatography with high concentration of organic solvent in mobile phase) concentration of organic solvent in mobile phase)	Ensure mobile phase compatibility with buffer concentration; decrease ionic strength and water-organic solvent ratio; premix mobile phase
	When injector disconnected from column - blockage in injector	Clean injector or replace rotor
Increasing Pressure	Blocked flow lines	Systematically disconnect components from detector end to column end to find blockage; replace or clean blocked component
	Particulate buildup at head of column	"Filter sample; use .5 micrometer in-line filter; disconnect and backflush column; replace inlet frit"
	Water-organic solvent systems - buffer precipitation	"Ensure mobile phase compatibility with buffer concentration; decrease ionic strength or water organic solvent ratio"
Broad peaks	Analytes eluted early due to sample overload	Dilute sample 1:10 and reinject
	Detector-cell volume too large	Use smallest possible cell volume consistent with sensitivity needs; use detector with no heat exchanger in system
	Injection volume too large	Decrease solvent strength of injection solvent to focus solute; inject smaller volume
	Large extra column volume	Use low- or zero-dead-volume endfittings and connectors; use smallest possible diameter of connecting tubing (<0.10 in. i.d.); connect tubing with matched fittings
	Mobile-phase solvent viscosity too high	Increase column temperature; change to lower viscosity solvent
	Peak dispersion in injector valve	Decrease injector sample loop size; introduce air bubble in front and back of sample in loop
	Poor column efficiency	Use smaller-particle-diameter packing, lower-viscosity mobile phase, higher column temperature, or lower flow rate
	Retention time too long	Use gradient elution or stronger isocratic mobile phase
	Sampling rate of data system too low	Increase sampling frequency.
	Slow detector time constant	Adjust time constant to match peak width
	Some peaks broad - late elution of analytes retained from previous injection	Flush column with strong solvent at end of run; end gradient at higher solvent concentration
Ghost peaks	Contamination	Flush column to remove contaminant; use HPLC-grade solvent
	Elution of analytes retained from previous injection	Flush column with strong solvent at end of run; end gradient at higher solvent concentration
	Ion-pair chromatography - upset equilibrium	Prepare sample in mobile phase; reduce injection volume
	Oxidation of trifluoroacetic acid in peptide mapping	Prepare trifluoroacetic acid solutions fresh daily; use antioxidant
	Reversed-phase chromatography - contaminated water	Check suitability of water by running different amounts through column and measure peak height of interferences as function of enrichment time; clean water by running it through old reversed-phase column; use HPLC-grade water.
	Unknown interferences in sample	Use sample cleanup or prefractionation before injection.
Negative peaks	Refractive index detection - refractive index of solute less than that of mobile phase	Reverse polarity to make peak positive
	UV-absorbance detection - absorbance of solute less than that of mobile phase	Use mobile phase with lower UV absorbance; if recycling solvent, stop recycling when recycled solvent affects detection

CHROMATOGRAPHY

Problem	Possible cause	Solution
Peak Doubling	Blocked Frit	Replace or clean frit; install 0.5-um porosity in-line filter between pump and injector to eliminate mobile-phase contaminants or between injector and column to eliminate sample contaminants
	Coelution of interfering compound	Use sample cleanup or prefractionation; adjust selectivity by changing mobile or stationary phase
	"Coelution of interfering compound from previous injection"	Flush column with strong solvent at end of run; end gradient at higher solvent concentration
	Column overloaded	Use higher-capacity stationary phase; increase column diameter; decrease sample amount
	Column void or channeling	Replace column, or, if possible, open top endfritting and clean and fill void with glass beads or same column packing; repack column
	Injection solvent too strong	Use weaker injection solvent or stronger mobile phase
	Sample volume too large	Use injection volume equal to one-sixth of column volume when sample prepared in mobile phase for injection
	Unswept injector flow path	Replace injector rotor
Peak Fronting	Channeling in column	Replace or repack column
	Column overloaded	Use higher-capacity stationary phase; increase column diameter; decrease sample amount
Tailing Peaks	Basic solutes - silanol interactions	Use competing base such as triethylamine; use a stronger mobile phase; use base-deactivated silica-based reversed-phase column; use polymeric column
	Beginning of peak doubling	See peak doubling
	Chelating solutes - trace metals in base silica	Use high purity silica-based column with low trace-metal content; add EDTA or chelating compound to mobile phase; use polymeric column
	Silica-based column - degradation at high pH	Use polymeric, sterically protected, or high-coverage reversed-phase column; install silica gel saturator column between pump and injector
	Silica-based column - degradation at high temperature	Reduce temperature to less than 50 °C
	Silica-based column - silanol interactions	Decrease mobile-phase pH to suppress silanol ionization; increase buffer concentration; derivatize solute to change polar interactions
	Unswept dead volume	Minimize number of connections; ensure injector rotor seal is tight; ensure all compression fittings are correctly seated
	Void formation at head of column	Replace column, or, if possible, open top end fitting and clean and fill in void with glass beads or same column packing; rotate injection valve quickly; use injection valve with pressure bypass; avoid pressure shock
Spikes	Bubbles in mobile phase	Degas mobile phase; use back-pressure restrictor at detector outlet; ensure that all fittings are tight
	Column stored without caps	Store column tightly capped; flush reversed-phase columns with degassed methanol

Correspond with other brand columns

Column	Supelco	Kromasil	Agilent	GL
C18-WP	Discovery RP-Amide C16		ZORBAX Rx C18	Inertsil ODS-EP
C18	SUPELCOSIL LC-18 Discovery C18	Kromasil C-18	ZORBAX Eclipse XDB-C18	Inertsil ODS-2
C18-BIO	Discovery BIO Wide Pore C18	Kromasil 300A C-18	ZORBAX 300SB-C18	Inertsil WP300 c18
C8	DISCOVERY C8	Kromasil C-8	ZORBAX Eclipse XDB-C8	Inertsil C8
C4		Kromasil C4		Inertsil C4
Phenyl	SUPELCOSIL LC-DP	Kromasil Phenyl	ZORBAX Eclipse XDBPhenyl	
Silica	SUPELCOSIL LC-Si	Kromasil SIL	ZORBAX Silica	Inertsil Sil
NH2	SUPELCOSIL LC-NH2	Kromasil NH2	ZORBAX NH2	Inertsil NH2
CN	SUPELCOSIL LC-CN	Kromasil CN	ZORBAX Eclipse XDB-CN	Inertsil CN-3

Column	Merck	Waters	Thermo
C18-WP		SymmetryShield C18	
C18	Puorsher STAR RP-18 endcapped	Symmetry C18	Hypersil ODS C18
C18-BIO	Lichrospher wp 300 RP-18e	Symmetry 300	Hypersil 300A C18
C8	Puorsher STAR RP-8 endcapped	Symmetry C8	Hypersil C8
C4		Spherisorb® C4	Hypersil GOLD C4
Phenyl		Spherisorb® Phenyl	Hypersil Phenyl-2
Silica	Lichrospher si 100	Spherisorb® W(Silica)	Hypersil Silica
NH2	Puorsher STAR NH2	Spherisorb® NH2	Hypersil NH2
CN	Lichrospher CN	Spherisorb® CN	Hypersil CN (CPS-2)

The USP liquid phase column summary

USP is United States Pharmacopoeia, provides a number of indicators for HPLC column packing:

USP	Packing Description	Recommend HPLC columns
L1	Octadecyl silane chemically bonded to porous silica or ceramic μ particles, 1.5 to 10 μ m in diameter, or a monolithic rod	C18, C18-WP, C18-BIO
L2	Octadecyl silane chemically bonded to silica gel of a controlled surface porosity that has been bonded to a solid spherical core, 30 to 50 μ m in diameter	C18 Packing
L3	Porous silica microparticles, 5 to 10 μ m in diameter	Silica
L4	Silica gel of controlled surface porosity bonded to a solid spherical core, 30 to 50 μ m in diameter	Silicycle packing
L7	Octyl silane chemically bonded to totally porous microsilica particles, 3 to 10 μ m in diameter	C8
L8	An essentially monomolecular layer of aminopropyl-silane chemically bonded to totally porous silica gel support, 10 μ m in diameter	NH2
L9	Irregular or spherical, totally porous silica gel having a chemically bonded, strongly acidic ation-exchange coating, 3 to 10 μ m in diameter	SCX
L10	Nitrile groups chemically bonded to porous silica microparticles, 3 to 10 μ m in diameter	CN, Shodex Silica 5CN
L11	Phenyl groups chemically bonded to porous silica microparticles, 3 to 10 μ m in diameter	Phenyl, Shodex Silica 5NPE
L12	Strong anion-exchange packing made by chemically bonding a quaternary amine to a solid silica spherical core, 30 to 50 μ m in diameter"	
L13	Trimethylsilane chemically bonded to porous silica microparticles, 3 to 10 μ m in diameter	Shodex Silica 5TMS
L14	Silica gel having a chemically bonded, strongly basic quaternary ammonium anion-exchange coating, 5 to 10 μ m in diameter.	SAX
L15	Hexyl silane chemically bonded to totally porous silica particles, 3 to 10 μ m in diameter	Spherisorb S5 C6
L16	Dimethyl silane chemically bonded to totally porous silica particles, 5 to 10 μ m in diameter	
L17	Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the hydrogen form, 7 to 11 μ m in diameter	Sep H-L, H-M, H-H
L18	Dimethyl silane chemically bonded to totally porous silica particles, 5 to 10 μ m in diameter	
L19	Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the calcium form, 9 μ m in diameter.	Sep Ca-L, Ca-M, Ca-H
L20	Dihydroxypropane groups chemically bonded to porous silica particles, 3 to 10 μ m in diameter.	Shodex PROTEIN KW-800
L21	A rigid, spherical styrene-divinylbenzene copolymer, 5 to 10 μ m in diameter.	"Transgenomic PRX-1, Shodex GPC KF-800,K-800, KD-800"
L22	A cation exchange resin made of porous polystyrene gel with sulfonic acid groups, about 10 μ m in size	Shodex ICY-521, SUGAR KS-800 series
L23	An ion exchange resin made of porous polymethacrylate or polyacrylate gel with quaternary ammonium groups, about 10 μ m in size	Shodex IEC QA-825
L24	A semi-rigid hydrophilic gel consisting of vinyl polymers with numerous hydroxyl groups on the matrix surface, 32 to 63 μ m in diameter	
L25	Packing having the capacity to separate compounds with a MW range from 100 to 5000 daltons (as determined by polyethylene oxide), applied to neutral, anionic, and cationic water- soluble polymers. A polymethacrylate resin base,crosslinked with poly-hydroxylated ether (surface contained some residual carboxyl functional groups) was found suitable.	Shodex OHpak SB-802 HQ Shodex OHpak SB-802.5 HQ, SB402.5
L26	Butyl silane chemically bonded to totally porous silica particles, 5 to 10 μ m in diameter	C4
L27	Porous silica particles, 30 to 50 μ m in diameter	Silicycle packing
L28	A multifunctional support, which consists of a high purity, 100 , spherical silica substrate that has been bonded with anionic (amine) functionality in addition to a conventional reversed phase C8 functionality"	
L29	Gamma alumina, reversed phase, low carbon percentage by weight, alumina-based polybutadiene spherical particles, 5 μ m diameter with a pore diameter of 80"	
L30	Ethyl silane chemically bonded to a totally porous silica particle, 3 to 10 μ m in diameter	
L31	A strong anion-exchange resin-quaternary amine bonded on latex particles attached to a core of 8.5 μ m macroporous particles having a pore size of 2000 Å and consisting of ethylvinylbenzene cross-linked with 55 % divinyl benzene	
L32	A chiral ligand-exchange packing- L-proline copper complex covalently bonded to irregularly shaped silica particles, 5 to 10 μ m in diameter	
L33	Packing having the capacity to separate proteins of 4,000 to 400,000 daltons. It is spherical, silica-based and processed to provide pH stability	Shodex PROTEIN KW-800 series Shodex KW400 series
L34	Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the lead form, about 9 μ m in diameter	Sep Pb-L, Pb-M, Pb-H
L35	A zirconium-stabilized spherical silica packing with a hydrophilic (diol-type) molecular monolayer bonded phase having a pore size of 150Å.	Agilent Zorbax GF-250
L36	3,5-dinitrobenzoyl derivative of L-phenylglycine covalently bonded to 5 μ m aminopropyl silica	
L37	Packing having the capacity to separate proteins by molecular size over a range of 2,000 to 40,000 Da. It is a polymethacrylate gel	Shodex OHpak SB-803 HQ, SB403
L38	Methacrylate-based size-exclusion packing for water-soluble samples	Shodex OHpak SB-802HQ
L39	Hydrophilic polyhydroxymethacrylate gel of totally porous spherical resin	Shodex Ohpak SB-800HQ, Shodex Rspak DM-614
L40	Cellulose tri-3,5-dimethylphenylcarbamate coated porous silica particles, 5 μ to 20 μ m in diameter	Regis Cell ®
L41	Immobilized α 1-acid glycoprotein on spherical silica particles	
L42	Octylsilane and octadecylsilane groups chemically bonded to porous silica particles,5 μ m in diameter	
L43	Pentafluorophenyl groups chemically bonded to silica particles, 5 to 10 μ m in diameter(5-10 μ m)	Supelco Discovery HSF5
L44	A multifunctional support, which consists of a high purity, 60 , spherical silica substrate that has been bonded with a cationic exchanger, sulfonic acid functionality in addition to a conventional reversed phase C8 functionality.	
L45	Beta cyclodextrin bonded to porous silica particles, 5 to 10 μ m in diameter	Shodex ORpak CDBS-453

CHROMATOGRAPHY

USP	Packing Description	Recommence HPLC columns
L46	Polystyrene/divinylbenzene substrate agglomerated with quaternary amine functionalized latex beads, 10 µm in diameter.	
L47	High capacity anion-exchange microporous substrate, fully functionalized with a trimethylamine group, 8 µm in diameter.	
L48	Sulfonated, cross-linked polystyrene with an outer layer of submicron, porous, anion-exchange microbeads, 15 µm in diameter.	
L49	A reversed-phase packing made by coating a thin layer of polybutadiene on to spherical porous zirconia particles, 3 to 10 µm in diameter.	Discovery Zr-PBD
L50	Multifunction resin with reversed-phase retention and strong anion-exchange functionalities. The resin consists of ethylvinylbenzene, 55 % cross-linked with divinylbenzene copolymer, 3 to 15 µm in diameter, and a surface area of not less than 350 m ² /g, substrate is coated with quaternary ammonium functionalized latex particles consisting of styrene cross-linked with divinylbenzene.	
L51	Amylose tris-3,5-dimethylphenylcarbamate-coated, porous, spherical, silica particles, 5 to 10 µm in diameter.	®
L52	A strong cation exchange resin made of porous silica with sulfopropyl groups, 5 to 10 µm in diameter.	SCX
L53	Weak cation-exchange resin consisting of ethylvinylbenzene, 55 % cross-linked with divinylbenzene copolymer, 3 to 15 µm diameter. Substrate is surface grafted with carboxylic acid and/or phosphoric acid functionalized monomers. Capacity not less than 500 µm in diameter.	
L54	"A size exclusion medium made of covalent bonding of dextran to highly cross-linked porous agarose beads, about 13 µm in diameter."	
L55	A strong cation exchange resin made of porous silica coated with polybutadiene-maleic acid copolymer, about 5 µm in diameter.	
L56	Isopropyl silane chemically bonded to totally porous silica particles, 3 to 10 µm in diameter	
L57	A chiral-recognition protein, ovomucoid, chemically bonded to silica particles, about 5 µm in diameter, with a pore size of 120 angstroms.	
L58	Strong cation-exchange resin consisting of sulphonated cross-linked styrene-divinylbenzene copolymer in the sodium form, about 7 to 11µm diameter	Sep Na-L, Na-M, Na-H, Transgenomic Coregel 87N
L59	Packing having the capacity to separate proteins by molecular weight over the range of 10 to 500kDa. It is spherical(10µm), silica-based,and processed to provide hydrophilic characteristics and pH stability	Shodex PROTEIN KW-800 series, Shodex KW400 series"
L60	Spherical, porous silica gel, 3 to 10 µm in diameter, surface has been covalently modified with palmitamidopropyl groups and endcapped.	C18-WP
L61	Hydroxide-selective, strong anion-exchange resin consisting of a highly cross-linked core of 13 µm microporous particles, pore size less than 10 , and consisting of ethylvinylbenzene cross-linked with 55 % divinylbenzene with a latex coating composed of 85 nm diameter microbeads bonded with alkanol quaternary ammonium ions (6 %).	
L62	C30 silane bonded phase on a fully porous spherical silica, 3 to 15 µm in diameter.	C30
L63	Glycopeptide teicoplanin linked through multiple covalent bonds to a 100 Å units spherical silica	
L64	Strongly basic anion exchange resin consisting of 8% crosslinked styrene divinylbenzene copolymer with a quaternary ammonium group in the chloride form, 45 to 180 µm in diameter	
L65	Strongly acidic cation exchange resin consisting of 8% sulfonated crosslinked styrene divinylbenzene copolymer with a sulfonic acid group in the hydrogen form, 63 to 250 µm in diameter	
L66	A crown ether coated on a 5 µm particle size silica gel substrate. The active site is (S)-18-crown-6 ether	
L67	Porous vinyl alcohol copolymer with a C18 alkyl group attached to the hydroxyl group of the polymer, 2 to 10 µm in diameter	Shodex Asahipak ODP-40 Shodex ET-RP1
L68	Spherical, porous silica, 10µm or less in diameter, the surface of which has been covalently modified with alkyl amide groups and not endcapped	
L69	Ethylvinylbenzene/divinylbenzene substrate agglomerated with quaternary amine functionalized 130nm latex beads, about 6.5µm in diameter	
L70	Cellulose tris(phenyl carbamate) coated on 5µm silica	
L71	Arigid, spherical polymetacrylate, 4 to 6 µm in diameter	Shodex RSpak DE-613
L72	(S)-phenylglycine and 3,5-dinitroaniline urea linkage covalently bonded to silica	
L73	A rigid, spherical polydivinylbenzene particle, 5 to 10 µm in diameter	
L74	A strong anion-exchange resin consisting of a highly cross-linked core of 7-µm macroporous particles having a 100 Angstroms average pore size and consisting of ethylvinylbenzene cross-linked with 55% divinylbenzene and an anion-exchange layer grafted to the surface, which is functionalized with alkyl quaternary ammonium ions.	
L75	A chiral-recognition protein, bovine serum albumin (BSA), chemically bonded to silica particles, about 7 µm in diameter, with a pore size of 300 Angstroms.	

Pressure unit conversion table

1atm = 1.01325bar

UNIT	Pa	KPa	MPa	bar	kgf/cm ²	mmH ₂ O	mmHg	p.s.i
Pa	1	10 ⁻³	10 ⁻⁶	10 ⁻⁵	10.2×10 ⁻⁶	101.97×10 ⁻³	7.5×10 ⁻³	0.15×10 ⁻³
KPa	10 ³	1	10 ⁻³	10 ⁻²	10.2×10 ⁻³	101.97	7.5	0.15
MPa	10 ⁶	10 ³	1	10	10.2	101.97×10 ³	7.5×10 ³	0.15×10 ³
bar	10 ⁵	10 ²	10 ⁻¹	1	1.02	10.2×10 ³	750.06	14.5
kgf/cm ²	98066.5	98.07	98.07×10 ⁻³	0.98	1	10.000	735.56	14.22
mmH ₂ O	9.806	9.807×10 ⁻³	9.807×10 ⁻⁶	98.07×10 ⁻⁶	10 ⁻⁴	1	73.56×10 ⁻³	1.42×10 ⁻³
mmHg	133.32	133.32×10 ⁻³	133.32×10 ⁻⁶	1.33×10 ⁻³	1.36×10 ⁻³	13.6	1	19.34×10 ⁻³
p.s.i	6894.76	6.89	6.89×10 ⁻³	68.95×10 ⁻³	70.31×10 ⁻³	703.07	51.71	1


Pressure unit conversion table

Solvent ①②	UV wavelength nm ③	Refractive index ④	Boiling point °C	Viscosity (cp 25 °C)	Polarity	Solubility ⑤	Dielectric constant 20 °C
Isooctane (*)	210	1.389	99	0.47	0.1	0.01	1.94
N-heptane (*)	200	1.385	98	0.4	0.2	0.01	1.92
N-hexane (*)	190	1.372	69	0.3	0.1	0.01	1.88
N-pentane (**)	210	1.355	36	0.22	0	0.01	1.84
Cyclohexane	210	1.423	81	0.9	0.1	0.012	2.02
Cyclopentane (*)	210	1.404	49	0.42	0.2	0.004	1.97
Carbon tetrachloride	265	1.457	77	0.9	1.6	0.008	2.24
Toluene	285	1.494	110	0.55	2.4	0.046	2.4
Xylene	290	1.493	138	0.6	2.5	unknown	2.3
Chlorobenzene	unknown	1.521	132	0.75	2.7	unknown	5.6
Benzene	280	1.498	80	0.6	2.7	0.07	2.3
Dichloromethane (**)	245	1.421	40	0.41	3.1	1.6	8.9
N-butanol	210	1.397	118	2.98	3.9	7.81	17.5
N-propanol	210	1.385	97	2.27	4	Miscible	20.3
Tetrahydrofuran(*)	220	1.405	66	0.55	4	Miscible	7.4
Ethyl acetate (*)	256	1.37	77	0.43	4.4	8.7	6.4
Isopropanol	210	1.384	82	2.3	4.3	Miscible	18.3
Chloroform (*)	245	1.443	61	0.53	4.1	0.815	4.8
Acetone (*)	330	1.356	56	0.3	5.4	Miscible	21.4
Ethanol	210	1.359	78	1.08	4.3	Miscible	24.6
Acetic acid	230	1.37	118	1.26	6	Miscible	6.2
Acetonitrile	210	1.341	82	0.34	6.2	Miscible	37.5
Methanol (*)	210	1.326	65	0.54	6.6	Miscible	32.7
Glycol	unknown	1.431	197	19.9	6.9	Miscible	37.7
Water	268	1.338	100	1	10.2	Miscible	80

- ① (*) means a low viscosity (<0.5cp), boiling point appropriate in (> 45 °C)
 ② (**) means small viscosity, low boiling point solvent.
 ③ Means approximate cutoff wavelength, when lower than this value, solvent is opaque.

- ④ Refractive index when 25 °C .
 ⑤ Percentage by weight of water at 20 °C when dissolved in a solvent, this value is useful in the liquid - solid chromatography.

Solvent miscibility



Application Index

Melamine in Milk Powder (according to GB/T 22388-2008)	13	VitB6.....	43	Malachite green and crystal violet aquatic ..	53
Tricyclic antidepressants	14	Coenzyme Q	43	Nitroaniline	53
Purine alkaloid	14	Psoralen	43	Polycyclic aromatic hydrocarbons (PAHs) (HJ 478-2009).....	54
Oligonucleotide	15	Loganin.....	43	Tetracyclines	54
Anti-HIV drugs	15	Paeonol	44	Watery nonvolatile pesticides	54
Tricyclic antidepressants	16	Cefixime	44	Bisphenol A.....	54
SulfaNo	16	Clarithromycin	44	Leather's phthalates.....	55
Hydrolysis bovine serum albumin.....	17	Paeoniflorin	44	Parabens in cosmetics.....	55
O-phthalic monoester acid.....	17	Carbamazepine.....	45	Phthalate monoester.....	55
Fat-soluble vitamins	20	Acetylacetone	45	Bromopyrene-C8	55
Steroid	20	Domiphen bromide.....	45	Nucleoside -1	56
Carbohydrates	21	Methotrexate	45	Oligonucleotide	56
Tocopherol isomers	21	Water-soluble vitamins.....	46	Nucleoside	56
Melamine.....	24	Fat-soluble vitamins.....	46	Hydrolysis of bovine serum albumin.....	57
Protein separation.....	27	Vitamin B	46	Synthetic peptide	57
Sorbitol and Mannitol	28	Citrus red No. 2 in juice.....	46	Protein sample	57
Comparison of different columns for separation of protein samples.....	33	Carbohydrate -1	47		
Series protein molecular weight calibration curve.....	33	Carbohydrate -2	47		
Tricyclic antidepressants	37	Carbohydrate -3	47		
Beta-blockers	37	Isomaltooligosaccharide	47		
Cough and cold medicine ingredients	37	Melamine.....	48		
Procainamide	37	Melamine in Milk Powder.....	48		
Anticholinergics	38	Melamine in raw milk (according to GB / T 22400-2008)	48		
Non-steroidal anti-inflammatory drugs	38	Furosine	48		
Glycyrrhizin	38	Vanillin and ethyl vanillin.....	49		
Matrine	38	Tocopherol isomers	49		
Anti-HIV drugs	39	Tocopherol isomers	49		
Tricyclic antidepressants	39	Tocopherol	49		
Steroids -1	39	Benzoic acid, sorbic acid	50		
Steroids -2	39	Sorbitol and mannitol	50		
Doxepin Hydrochloride	40	Sudan in chili sauce.....	50		
Cefotaxime valerate.....	40	Synthetic colorants	50		
Deoxyschizandrin schisandrin B	40	Carbamate pesticide in pepper.....	51		
Calcium pantothenate	40	Fungicides	51		
Cefuroxime Sodium	41	Glyphosate	51		
Taurine	41	Clethodim	51		
Rifampicin and related substances	41	Quinolones	52		
Melatonin.....	41	Sulfa	52		
Day ephedra.....	42	Sulfa drugs in feed	52		
Berberine	42	Clenbuterol.....	52		
Aspirin C in propolis.....	42	Fluoroquinolones	53		
Nicotinamide	42	Norfloxacin	53		

Brief introduction

Currently, HPLC is widely used in the chemical, biological and pharmaceutical field. These series include silica and polymer matrix columns, both analysis and preparative columns, to meet needs of customers in various fields.

Silica-based analytical column

HPLC columns

The Columns base on high-purity silica gel, using unique bonding technique, with excellent peak shape, better selectivity, sensitivity and reproducibility. With low content of matrix metal, the columns show perfect peak shape for all types of analytes. Different types of bonded phases provide more flexibility for method development. To ensure excellent column performance and long column life, we comply with strict production process in manufacturing and have a strict quality control for HPLC columns.

- Suitable for all types of samples
- Excellent column reproducibility
- A variety of bonded phases

The packings information:

Packings	C18-WP	C18	C18-BIO	C8	C4	Phenyl	CN	Diol
Particle diameter (μm)	3 and 5	5 and 10	5	3 and 5	5	5	3 and 5	3 and 5
Pore size(Å)	100	120	300	120	300	120	120	120
Pore volume (mL/g)	1.1	1.0	0.9	1.0	0.9	1.0	1.0	1.0
Endcapped	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Specific surface area(m ² /g)	450	300	100	300	100	300	300	300
Metallic impurities (ppm)	<10	<10	<10	<10	<10	<10	<10	<10
Carbon content	17%	17%	8%	10%	3%	11%	7.5%	8.8%
pH range	1.5 - 10	2 - 8	1.5 - 11	2 - 8	2 - 8	2 - 8	2.5 - 8	2.5 - 8
Temperature range (°C)	20 - 60	20 - 60	20 - 60	20 - 60	20 - 60	20 - 60	20 - 60	20 - 60

Packings	NH ₂	Silica	SAX	SCX	HILIC	HILIC(2)	HILIC(3)	30
Particle diameter (μm)	3 and 5	3 and 5	3 and 5	3 and 5	3 and 5	3 and 5	5	3 and 5
Pore size(Å)	120	120	120	120	120	120	120	120
Pore volume (mL/g)	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0
Specific surface area(m ² /g)	300	300	300	300	300	300	300	450
Metallic impurities (ppm)	<10	<10	<10	<10	<10	<10	<10	<10
Carbon content	4%	0%	16%	11%	8.6%	8%	16%	20%
pH range	2 - 8	2 - 8	2 - 8	2 - 8	1.5 - 8	1.5 - 8	1.5 - 8	2 - 8
Temperature range (°C)	20 - 60	20 - 60	20 - 60	20 - 60	20 - 60	20 - 60	20 - 60	20 - 60




C18-WP

[Recommended for Method Development, fit for a variety of mobile phase conditions]

C18-WP use high purity of spherical silica matrix and have excellent stability. C18-WP can use 100% pure water as mobile phase for separation of acidic, neutral and basic organic compound, as well as many drugs and peptides etc. A variety of specifications, from analytical to preparative scale can be provided.

- bonded C18 groups
- pH stability range: 1.5-10
- Suitable for 100% water mobile phase
- Strong retain for polar substances
- Symmetrical peak shape for Alkaline substances
- High specific surface area, suitable for high load



PH stability


Stability of low pH

In the low pH mobile phase, the main reason for short column life is drop of chemical bonded group from silica gel by hydrolysis. Hydrolysis leads to changing retention time of the analyte, short lifetime and poor reproducibility.

The following figure shows C18-WP stability under the conditions of pH 1.5 mobile phase.

Low pH tolerance (pH 1.5)

Column	C18-WP, 4.6 x 150 mm, 5µm
Mobile phase	Acetonitrile: 0.1% trifluoroacetic acid (pH 1.5) (50/50)
Flow rate	1.0 mL / min
Detection	UV 254 nm
Column temperature	30 °C
Sample	toluene




Stability of high pH

In the high pH mobile phase, silica matrix is gradually dissolved. Ordinary pH range of silica-based columns is 2-8. When the pH of mobile phase is more than 8, silica gel is dissolved speedily, and column life is very short. C18-WP columns can protect silica matrix to have a longer life in high pH conditions, due to unique bonding and endcapped technology.

High pH tolerance (pH 11.0)

Column	C18-WP, 4.6 x 150 mm, 5µm
Mobile phase	Methanol: 0.5% aqueous ammonia (pH 11.0) (20/80)
Flow rate	1.0 mL / min
Detection	UV 254 nm
Column temperature	30 °C
Sample	Phthalatedipropyl



100% stability of the aqueous phase

Usually silica-based reversed-phase column can not be used in high proportion of water mobile phase conditions, and organic phase in the mobile phase must be maintained more than 5%. This may limit some polar compounds' separation in reversed-phase conditions. The reason is hydrophobic collapse.

"Hydrophobic collapse" is a phenomenon that reversed-phase column loss the ability of retaining compounds in a mobile phase with a very high water content. Due to the hydrophobic interaction of functional groups, the surface of the stationary phase cannot be wet by the mobile phase and hydrophobic chains fold up.


According to the research, the hydrophobic collapse generally occurs when restarting of the mobile phase after stopping pump. The experiment can verify whether a column is compatible with pure water. Test column efficiency at first, and wash the column with 100% water mobile phase at 1.0 mL/min for 2h. Then slow down the flow rate to zero and stop pump for 1h. Columns are washed with 100% water mobile phase again and tested for column efficiency the second time.

Compare the difference of retention before and after stopping pump.


Test Condition

Column	C18-WP, 4.6 x 150 mm, 5µm
Mobile phase	methanol: water (70/30)
Flow rate	1.0 mL / min
Detection	UV 254 nm
Column temperature	30 ° C
Sample	1. Uracil 2. Toluene 3. Naphthalene

(A)



(B)



100% water mobile phase compatibility


A before pump stopping

B after pump stopping

Melamine in Milk Powder (according to GB/T22388-2008) No.03215

1. Impurities in Milk

2. Melamine



Column: C18-WP 4.6 x 150mm, 5µm (HCA050U046X15072A)
 Mobile phase: 10 mM hexane sulfonate +10 mM citric acid buffer solution / acetonitrile (90/10)
 Flow rate: 1.0 mL/min
 Detection: 240 nm
 Column temperature: 40 °C

C18-WP


Ordering Information:

Product Code	Particle size	diameter x length
HCA030U021X05073A	3µm	2.1 x 50mm
HCA030U021X10073A	3µm	2.1 x 100mm
HCA030U021X15073A	3µm	2.1 x 150mm
HCA030U021X20073A	3µm	2.1 x 200mm
HCA030U021X25073A	3µm	2.1 x 250mm
HCA030U046X05073A	3µm	4.6 x 50mm
HCA030U046X10073A	3µm	4.6 x 100mm
HCA030U046X15073A	3µm	4.6 x 150mm
HCA030U046X20073A	3µm	4.6 x 200mm
HCA030U046X25073A	3µm	4.6 x 250mm
HCA050U021X05072A	5µm	2.1 x 50mm
HCA050U021X10072A	5µm	2.1 x 100mm
HCA050U021X15072A	5µm	2.1 x 150mm
HCA050U021X20072A	5µm	2.1 x 200mm
HCA050U021X25072A	5µm	2.1 x 250mm
HCA050U046X05072A	5µm	4.6 x 50mm
HCA050U046X10072A	5µm	4.6 x 100mm
HCA050U046X15072A	5µm	4.6 x 150mm
HCA050U046X20072A	5µm	4.6 x 200mm
HCA050U046X25072A	5µm	4.6 x 250mm

C18

[Conventional C18 column]

- Bonded C18 groups
- High-purity silica, metal content <10ppm
- Less hydrophobic than C18-WP, with different selectivity
- Economic column




Based on high purity spherical silica, C18 column is good at separating a variety of compounds. It is a typical economic column with high price ratio, as well as long column lifetime. For most analytes, the retention times are shorter than that of C18-WP columns of same specifications.

Purine alkaloid

No.03217

- | | |
|----------------|------------|
| 1.Theobromine | 3.Caffeine |
| 2.Theophylline | 4.Phenol |




Column: C18 4.6 × 150mm, 5μm (HCA050U046X15071A)
 Mobile phase: methanol / water (25/75)
 Flow rate: 1.0 mL/min
 Detection: 254 nm
 Column temperature: 40 °C

Tricyclic antidepressants

No.03216

- 1.Protriptyline
- 2.Nortriptyline
- 3.Amitriptyline



Column: C18 4.6 × 150mm, 5μm (HCA050U046X15071A)
 Mobile phase: methanol / 20 mM K₂HPO₄ buffer (pH 7.0) (80/20)
 Flow rate: 1.0 mL/min
 Detection: 254 nm
 Column temperature: 40 °C

Ordering Information:


Particle size diameter × length Item No.

Product Code	Particle size	diameter × length
HCA050U021X05071A	5μm	2.1 × 50mm
HCA050U021X10071A	5μm	2.1 × 100mm
HCA050U021X15071A	5μm	2.1 × 150mm
HCA050U021X20071A	5μm	2.1 × 200mm
HCA050U021X25071A	5μm	2.1 × 250mm
HCA050U046X05071A	5μm	4.6 × 50mm
HCA050U046X10071A	5μm	4.6 × 100mm
HCA050U046X15071A	5μm	4.6 × 150mm
HCA050U046X20071A	5μm	4.6 × 200mm
HCA050U046X25071A	5μm	4.6 × 250mm
HCA100U046X15074A	10μm	4.6 × 150mm
HCA100U046X25074A	10μm	4.6 × 250mm

C18-BIO

[Applicable to macromolecules]

- Bonded C18 groups
- 300Å pore size, fit for macromolecules separation, such as peptides
- High column efficiency and long lifetime
- Stable in the range of pH 1.5-11




300Å pore size, highpurity silica, high density bonding, and completely endcapped, make C18-BIO able to separatelarge moleculars, especially proteins and polypeptides.

Applications:

Oligonucleotide

No.03218

- | | |
|--------------|--------------|
| 1.CAAGACGCAA | 4.CCCTGAACAA |
| 2.CAACCAACGT | 5.CGTGTATTGG |
| 3.GGTGATCAAC | 6.GGTCCTATAC |




Column: C18-BIO 4.6 × 150mm, 5µm (HCA050U046X15078A)
 Mobile phase: A: 50 mM NaH₂PO₄ buffer solution (pH 7.0); B: acetonitrile 0min B: 5%; 20min B: 15%
 Flow rate: 1.0 mL/min
 Detection: 260 nm
 Column temperature: 25 °C

Anti-HIV drugs

No.03217

- | | |
|---------------|-------------------|
| 1.Theobromine | 3.AZT-Glucuronide |
| 2.d4T | 4.AZT |



Min

Column: C18-BIO 4.6 × 150mm, 5µm (HCA050U046X15078A)
 Mobile phase: methanol / 20 mM NH₄H₂PO₄ buffer (10/90)
 Flow rate: 1.0 mL/min
 Detection: 260 nm
 Column temperature: 35 °C


Ordering Information:

Product Code	Particle size	diameter × length
HCA050U021X05078A	5µm	2.1 × 50mm
HCA050U021X10078A	5µm	2.1 × 100mm
HCA050U021X15078A	5µm	2.1 × 150mm
HCA050U021X20078A	5µm	2.1 × 200mm
HCA050U021X25078A	5µm	2.1 × 250mm
HCA050U046X05078A	5µm	4.6 × 50mm
HCA050U046X10078A	5µm	4.6 × 100mm
HCA050U046X15078A	5µm	4.6 × 150mm
HCA050U046X20078A	5µm	4.6 × 200mm
HCA050U046X25078A	5µm	4.6 × 250mm

C8

[High resolution, rapid analysis]

- Bonded C8 groups
- Better resolution than C18 group for medium polarity subjects, and short retention time for non-polar compounds
- Good peak shapes for acidic, basic, and neutral substances
- Long column life and good repeatability




C8 offers less degree of hydrophobic selectivity compared to C18. C8 is a better choice if need to save time and achieve rapid analysis in the same chromatographic condition on octadecyl bonded phase.

Applications

SulfaNo

No.03221

- 1.Sulfonamide
- 2.Sulfisomidine
- 3.Sulfadiazine
- 4.Sulfamethazine




Column: C8 4.6 × 150mm, 5μm (HCA050U046X15075A)
 Mobile phase: acetonitrile / 0.1% H₃PO₄ buffer (10/90)
 Flow rate: 1.0 mL/min
 Detection: 254 nm
 Column temperature: 40 °C

Tricyclic antidepressants

No.03220

- 1.Protriptyline
- 2.Nortriptyline
- 3.Amitriptyline



Column: C8 4.6 × 150mm, 5μm (HCA050U046X15075A)
 Mobile phase: methanol / 20 mM K₂HPO₄ buffer (pH 7.0) (80/20)
 Flow rate: 1.0 mL/min
 Detection: 254 nm
 Column temperature: 40 °C


Ordering Information:

Product Code	Particle size	diameter × length
HCA030U021X05065A	3μm	2.1 × 50mm
HCA030U021X10065A	3μm	2.1 × 100mm
HCA030U021X15065A	3μm	2.1 × 150mm
HCA030U021X20065A	3μm	2.1 × 200mm
HCA030U021X25065A	3μm	2.1 × 250mm
HCA030U046X05065A	3μm	4.6 × 50mm
HCA030U046X10065A	3μm	4.6 × 100mm
HCA030U046X15065A	3μm	4.6 × 150mm
HCA030U046X20065A	3μm	4.6 × 200mm
HCA030U046X25065A	3μm	4.6 × 250mm
HCA050U021X05075A	5μm	2.1 × 50mm
HCA050U021X10075A	5μm	2.1 × 100mm
HCA050U021X15075A	5μm	2.1 × 150mm
HCA050U021X20075A	5μm	2.1 × 200mm
HCA050U021X25075A	5μm	2.1 × 250mm
HCA050U046X05075A	5μm	4.6 × 50mm
HCA050U046X10075A	5μm	4.6 × 100mm
HCA050U046X15075A	5μm	4.6 × 150mm
HCA050U046X20075A	5μm	4.6 × 200mm
HCA050U046X25075A	5μm	4.6 × 250mm

C4

[Low hydrophobic reverse phase, rapid analysis]

- Bonded C4 groups
- 300Å pore size, fit for macromolecules separation
- Rapid analysis
- High column efficiency and excellent peak shape




Retention times are shorter than on C8 and C18 phases. 300Å pore size is suitable for analysis of biological samples.

Phenyl

[Analysis for compounds with cyclic structure]


- Bonded phenylpropyl groups
- Interactions with π - π of aromatic compound
- Unique selectivity for compounds with cyclic structure
- good reproducibility



Phenyl column bonded phenylpropyl group, with surface coverage is 3.0 $\mu\text{mol}/\text{m}^2$. Phenyl exhibits a unique selectivity for aromatic compounds, due to a possibility for π - π interactions between the phenyl bonded phase and the solute.

Hydrolysis bovine serum albumin

No.03222




Column: C4 4.6 × 250mm, 5 μm (HCA050U046X15079A)
 Mobile phase: A: 0.09% TFA; B: 0.085% TFA + 80% acetonitrile
 0min B 5%; 5min B 5%; 35min B 50%; 45min B 100%
 Flow rate: 1.0 mL/min
 Detection: 214 nm
 Column temperature: 25 °C

Ordering Information:

Product Code	Particle size	diameter × length
HCA050U021X05079A	5 μm	2.1 × 50mm
HCA050U021X10079A	5 μm	2.1 × 100mm
HCA050U021X15079A	5 μm	2.1 × 150mm
HCA050U021X20079A	5 μm	2.1 × 200mm
HCA050U021X25079A	5 μm	2.1 × 250mm
HCA050U046X05079A	5 μm	4.6 × 50mm
HCA050U046X10079A	5 μm	4.6 × 100mm
HCA050U046X15079A	5 μm	4.6 × 150mm
HCA050U046X20079A	5 μm	4.6 × 200mm
HCA050U046X25079A	5 μm	4.6 × 250mm

O-phthalic monoester acid

No.03223



Column: Phenyl 4.6 × 150mm, 5 μm (HCA050U046X15037A)
 Mobile phase: acetonitrile / water / acetic acid(45/55/0.2)
 Flow rate: 0.8 mL/min
 Detection: 228 nm
 Column temperature: 25 °C

Ordering Information:

Product Code	Particle size	diameter × length
HCA050U021X05037A	5 μm	2.1 × 50mm
HCA050U021X10037A	5 μm	2.1 × 100mm
HCA050U021X15037A	5 μm	2.1 × 150mm
HCA050U021X25037A	5 μm	2.1 × 250mm
HCA050U046X05037A	5 μm	4.6 × 50mm
HCA050U046X10037A	5 μm	4.6 × 100mm
HCA050U046X15037A	5 μm	4.6 × 150mm
HCA050U046X25037A	5 μm	4.6 × 250mm

C30


[Applicable to Carotenoid Separation]

- Unique C30 bonded phase, offering diverse selectivity
- High shape selectivity for structurally similar isomers

C30 is bonded with unique C30 functional groups, suitable for the separation of polar substances (such as sugars and nucleic acids) and lipophilic compounds (such as vitamin E and carotenoids).


Determination of vitamins A and E

GB5009.82—2016



Column: C30, 4.6 x 250mm, 5μm (HCA050U046X25052A)
 Flow rate: 0.8ml/min
 Column temperature: 20°C
 Injection volume: 10μl
 Detector: UV, vitamin A at 325nm; vitamin E at 294nm
 Mobile phase: A: Methanol B: Water
 Gradient conditions:

Time (min)	Flow Rate	Water(%)	Methanol(%)
0.0	0.8	4	96
12.0	0.8	4	96
12.5	1.0	4	96
17.0	1.0	0	100
30.0	1.0	0	100
31.0	0.8	4	96
33.0	0.8	4	96



Column: C30, 4.6 x 250mm, 3μm (HCA030U046X25053A)
 Flow rate: 0.8ml/min
 Column temperature: 20°C
 Injection volume: 10μl
 Detector: UV, 294nm
 Mobile phase: 100% Methanol


Determination of Lutein

GB5009.248—2016


Column: C30, 4.6 x 250mm, 5μm (HCA050U046X25052A)
 Flow rate: 1.0ml/min
 Column temperature: 30°C
 Injection volume: 50μl
 Detector: UV, vitamin A at 325nm; vitamin E at 294nm
 Mobile phase: A:Methanol : Water (88:12, volume ratio, containing 0.1% BHT)
 B:Tert-butyl methyl ether(containing 0.1% BHT)
 Gradient conditions

Time (min)	A (%)	B (%)
0	85	15
5	85	15
18	60	40
18.1	85	15
25	85	15


Standard Curve



Sample Chromatogram:



Sample Chromatogram:




Ordering Information:

Packings	Product Code	diameter × length
C30	HCA050U046X25052A	4.6 × 250mm,5μm
C30	HCA030U046X25053A	4.6 × 250mm,5μm
C30 Guard Cartridge Kit	HCA050U040X02052KA	1 Holder and 1 Cartridge,5 μm,4.0 *20 mm
C30 Guard Cartridge	HCA050U040X02052A	4.0×20mm, 5μm

C18 HPLC Column

[Applicable to Carotenoid Separation]

- High column efficiency, low column pressure
- High-throughput rapid analysis
- Compatible with UHPLC and HPLC
- High reproducibility




Shell is a nucleoparticle silica chromatography column with a particle size of 2.6 μm . It consists of a solid spherical core with a diameter of 1.6 μm and a porous shell layer with a thickness of 0.5 μm . It possesses equal column efficiency and separation performance to sub-2 μm UHPLC columns while only requiring half the column pressure. This allows it to be used on both UHPLC and HPLC instruments, enabling ultra-fast separation with shorter analysis times and higher separation efficiency.solute.

HTraditional Chinese Medicine Pesticide Residues - Pharmacopoeia

- | | |
|---------------|---------------|
| 1.CAAGACGCAA | 4.CCCTGAACAA |
| 2.CAACCAACGT | 5. CGTGTATTGG |
| 3. GGTGATCAAC | 6. GGTCCATATC |

- | | | |
|------------------------------|-------------------------------------|-------------------------------------|
| 1.Methamidophos | 11.Benxlithium sulfoxide | 21.Benxlithium |
| 2.Temidiphos ethyl sulfoxide | 12.Metsulfuron methyl | 22.Fenthion |
| 3.Carbaryl | 13.Carbofuran | 23.Terbutyl thiophosphate sulfoxide |
| 4.Temidiphospropyl sulfoxide | 14.Chlorsulfuron | 24.Clorophene |
| 5.Ethyl parathion | 15.Metyl para-phosphorus sulfoxide | 25.Thioclorthophen |
| 6.3-Hydroxycarbofuran | 16.Aminopropyl sulfuron | 26.Methyl isothiophosphate |
| 7.Thiophos | 17.Ethyl parathion | 27.Ciclophen |
| 8.Benxlithium thiosulfate | 18.Terbutyl thiophosphate sulfoxide | 28.Dieldrin |
| 9.Phosphanil | 19.Metyl para-phosphorus sulfoxide | 29.Fenvalerate |
| 10.Temidiphos | 20.Nitrifurazone | 30.Metyl para-phosphorus |




Column: C18 2.1 \times 100mm, 2.6 μm (HCA026U021X1001A)
 A phase 0.1% formic acid solution (containing 5mmol/L ammonium formate)
 Mobile phase B phase acetonitrile - 0.1% formic acid solution (containing 5mmol/L ammonium formate) (95:5) wacetic acid(45/55/0.2)
 Flow rate: 0.3ml/min
 Temperature 40°C
 Ionization mod ESI, positive ion mode analysis
 modetemperature: Multi-reaction monitoring (MRM)

Time(min)	A(%)	B(%)
0~1	70	30
1~12	70~0	30~100
12~14	0	100

HTraditional Chinese Medicine Pesticide Residues - Pharmacopoeia

- | | |
|-----------------|--------------------------|
| 1.Sulfadiazine | 4.Sulfamethoxypyridazine |
| 2.Sulfathiazole | 5.Sulfamethoxazole |
| 3.Sulfamerazine | 6.Sulfaquinoxaline |



Column: C18, 2.1mm \times 50mm, 2.6 μm , (HCA026U021X05011A)
 Mobile phase A: Acetonitrile ; B: 0.1% Formic Acid
 Flow rate: 1.0ml/min
 Detector: 254nm
 Temperature 40°C

Time(min)	A(%)	B(%)
0	10	90
2	90	10

Ordering Information:

Packings	Product Code	Particle size	diameter \times length
C18	HCA026U021X05011A	2.6 μm	2.1 \times 50mm
C18	HCA026U021X10011A	2.6 μm	2.1 \times 100mm
C18	HCA026U046X10011A	2.6 μm	4.6 \times 100mm
C18	HCA026U046X05011A	2.6 μm	4.6 \times 50mm

More Information:

Packings	Product Code	Particle size	diameter \times length
PFP	HCA026U021X10013A	2.6 μm	4.6 \times 50mm

Silica

[Non- bonded silica, normal phase]


- Spherical silica, non- bonded
- For non-polar and medium polar organic compounds
- Ultra pure, low metal impurity
- Symmetrical peak shape

None bonded high-purity silica, metal impurity content <10ppm, high mechanical strength. Silica is fit for separation of non-polar and media polar organic compounds to achieve sharp peak shape and high reproducibility for columns.

Fat-soluble vitamins

No.03224

- 1.Vitamin A palmitate
- 2.Vitamin K1
- 3.Vitamin E
- 4.Vitamin K3
- 5.Vitamin D




Column: Silica 4.6 × 150mm, 5μm (HCA050U046X15076A)
 Mobile phase: n-hexane / chloroform (60/40)
 Flow rate: 1.0 mL/min
 Detection: 254 nm
 Column temperature: 25 °C

Steroid

No.03225

- 1.Estrone
- 2.Estradiol
- 3.Estriol



Column: Silica 4.6 × 150mm, 5μm (HCA050U046X15076A)
 Mobile phase: n-hexane / ethanol (85/15)
 Flow rate: 1.0 mL/min
 Detection: 270 nm
 Column temperature: 40 °C


Ordering Information:

Product Code	Particle size	diameter × length
HCA030U021X05066A	3μm	2.1 × 50mm
HCA030U021X10066A	3μm	2.1 × 100mm
HCA030U021X15066A	3μm	2.1 × 150mm
HCA030U021X20066A	3μm	2.1 × 200mm
HCA030U021X25066A	3μm	2.1 × 250mm
HCA030U046X05066A	3μm	4.6 × 50mm
HCA030U046X10066A	3μm	4.6 × 100mm
HCA030U046X15066A	3μm	4.6 × 150mm
HCA030U046X20066A	3μm	4.6 × 200mm
HCA030U046X25066A	3μm	4.6 × 250mm
HCA050U021X05076A	5μm	2.1 × 50mm
HCA050U021X10076A	5μm	2.1 × 100mm
HCA050U021X15076A	5μm	2.1 × 150mm
HCA050U021X20076A	5μm	2.1 × 200mm
HCA050U021X25076A	5μm	2.1 × 250mm
HCA050U046X05076A	5μm	4.6 × 50mm
HCA050U046X10076A	5μm	4.6 × 100mm
HCA050U046X15076A	5μm	4.6 × 150mm
HCA050U046X20076A	5μm	4.6 × 200mm
HCA050U046X25076A	5μm	4.6 × 250mm

NH₂

[Both Normal and reverse phase mode]

- Bonded aminopropyl group
- Suitable for normal and reverse phase mode
- Separate sugars in reverse mode




Aminopropyl stationary phase serves as a weak anion exchanger and offer polar selectivity under reversed phase and normal phase conditions.

Tocopherol isomers

No.03227

1. α -Tocopherol
2. α -Tocopherol
3. α -Tocopherol
4. α -Tocopherol




Min

Column: NH₂ 4.6 x 150mm, 5μm (HCA050U046X15077A)
 Mobile phase: n-hexane / ethyl acetate (70/30)
 Flow rate: 1.0 mL/min
 Detection: 295 nm
 Column temperature: 40 °C

Carbonhydrates

No.03226

1. Glucose
2. Maltose
3. Maltotriose
4. Maltotetraose
5. Matlopentaos



Min

Column: NH₂ 4.6 x 150mm, 5μm (HCA050U046X15077A)
 Mobile phase: acetonitrile / water (50/50)
 Flow rate: 1.0 mL/min
 Detection: RID
 Column temperature: 40 °C


Ordering Information:

Product Code	Particle size	diameter x length
HCA030U021X05067A	3μm	2.1 x 50mm
HCA030U021X10067A	3μm	2.1 x 100mm
HCA050U021X15067A	3μm	2.1 x 150mm
HCA030U021X20067A	3μm	2.1 x 200mm
HCA030U021X25067A	3μm	2.1 x 250mm
HCA030U046X05067A	3μm	4.6 x 50mm
HCA030U046X10067A	3μm	4.6 x 100mm
HCA030U046X15067A	3μm	4.6 x 150mm
HCA030U046X20067A	3μm	4.6 x 200mm
HCA030U046X25067A	3μm	4.6 x 250mm
HCA050U021X05077A	5μm	2.1 x 50mm
HCA050U021X10077A	5μm	2.1 x 100mm
HCA050U021X15077A	5μm	2.1 x 150mm
HCA050U021X20077A	5μm	2.1 x 200mm
HCA050U021X25077A	5μm	2.1 x 250mm
HCA050U046X05077A	5μm	4.6 x 50mm
HCA050U046X10077A	5μm	4.6 x 100mm
HCA050U046X15077A	5μm	4.6 x 150mm
HCA050U046X20077A	5μm	4.6 x 200mm
HCA050U046X25077A	5μm	4.6 x 250mm

CN

[Can be used for normal or reverse phase separation]

- Bonded cyanopropyl
- Can be used for normal or reverse phase separation
- High column efficiency and good reproducibility



CN is cyanide propyl bonded silica column with n-electron interaction and unshared electron pair hydrogen bonding. Can be used for both reverse phase and normal phase mode. When used for reverse mode, having different selectivity from C18 and C8 columns; when used for normal phase mode, retention is lower retention than non-bonded silica gel column.


Ordering Information:

Product Code	Particle size	diameter × length
HCA030U021X15034A	3μm	2.1 × 150mm
HCA050U021X05033A	5μm	2.1 × 50mm
HCA050U021X10033A	5μm	2.1 × 100mm
HCA050U021X15033A	5μm	2.1 × 150mm
HCA050U021X20033A	5μm	2.1 × 200mm
HCA050U021X25033A	5μm	2.1 × 250mm
HCA050U046X05033A	5μm	4.6 × 50mm
HCA050U046X10033A	5μm	4.6 × 100mm
HCA050U046X15033A	5μm	4.6 × 150mm
HCA050U046X20033A	5μm	4.6 × 200mm
HCA050U046X25033A	5μm	4.6 × 250mm

Diol

[Suitable for normal phase separation]

- Bonded group of 1,2 - dihydroxy-propyl ether propionate
- Normal phase separation
- Good reproducibility



Diol bonded 1,2 - dihydroxy-propyl ether propionate group, coverage of 4.0 micromol / m², can interact with polar compounds. diol is able to distinguish compounds from slight difference, also can separate biological molecules based on size exclusion mechanism.


Ordering Information:

Product Code	Particle size	diameter × length
HCA030U021X05036A	3μm	2.1 × 50mm
HCA030U021X10036A	3μm	2.1 × 100mm
HCA030U021X15036A	3μm	2.1 × 150mm
HCA030U021X25036A	3μm	2.1 × 250mm
HCA030U046X05036A	3μm	4.6 × 50mm
HCA030U046X10036A	3μm	4.6 × 100mm
HCA030U046X15036A	3μm	4.6 × 150mm
HCA030U046X25036A	3μm	4.6 × 250mm
HCA050U021X05035A	5μm	2.1 × 50mm
HCA050U021X10035A	5μm	2.1 × 100mm
HCA050U021X15035A	5μm	2.1 × 150mm
HCA050U021X25035A	5μm	2.1 × 250mm
HCA050U046X05035A	5μm	4.6 × 50mm
HCA050U046X10035A	5μm	4.6 × 100mm
HCA050U046X15035A	5μm	4.6 × 150mm
HCA050U046X25035A	5μm	4.6 × 250mm

SAX

[Suitable for analysis of acidic substances]

- Strong anion exchange mode
- Suitable for analysis of acidic substances, including nucleotide and organic acids etc.
- High column efficiency, high batch stability, good columns reproducibility
- To adjust retention time of the analytes by changing buffer concentration of mobile phase
- Stable in high proportion of water mobile phase



SAX column are bonded quaternary ammonium strong anion- exchange group in the high-purity silica matrix, having mixed chemical structure of quaternary ammonium and phenyl functional groups. This mixed-mode by strong anion exchange phase and hydrophobic phase is suitable for separation of aromatic or aliphatic carboxylic acids, sulfonic acids, nucleotides and acids etc.


Ordering Information:

Product Code	Particle size	diameter × length
HCA030U021X05020A	3μm	2.1 × 50mm
HCA030U021X10020A	3μm	2.1 × 100mm
HCA030U021X15020A	3μm	2.1 × 150mm
HCA030U021X25020A	3μm	2.1 × 250mm
HCA030U046X05020A	3μm	4.6 × 50mm
HCA030U046X10020A	3μm	4.6 × 100mm
HCA030U046X15020A	3μm	4.6 × 150mm
HCA030U046X25020A	3μm	4.6 × 250mm
HCA050U021X05021A	5μm	2.1 × 50mm
HCA050U021X10021A	5μm	2.1 × 100mm
HCA050U021X15021A	5μm	2.1 × 150mm
HCA050U021X25021A	5μm	2.1 × 250mm
HCA050U046X05021A	5μm	4.6 × 50mm
HCA050U046X10021A	5μm	4.6 × 100mm
HCA050U046X15021A	5μm	4.6 × 150mm
HCA050U046X25021A	5μm	4.6 × 250mm

SCX

[Suitable for analysis of alkaline substances]

- Strong cation exchange mode
- Fit for analysis of alkaline substances, especially amines
- high column efficiency, stable batch, good columns reproducibility



SCX is benzenesulfonic acid bonded silica, having mixed chemical structure of sulfonic acid group and phenyl group. SCX is mixed mode of strong cation exchange phase and hydrophobic phase. Not only can be used for separation of cationic / basic and nitrogenous compounds, but also give appropriate reservation for a variety of weak cation, neutral organic compound. SCX is used for separation and determination of amines and polyamine compounds, such as alkaloids, peptides and components in cold medicines.

Ordering Information:

Product Code	Particle size	diameter × length
HCA030U021X05022A	3μm	2.1 × 50mm
HCA030U021X10022A	3μm	2.1 × 100mm
HCA030U021X15022A	3μm	2.1 × 150mm
HCA030U021X25022A	3μm	2.1 × 250mm
HCA030U046X05022A	3μm	4.6 × 50mm
HCA030U046X10022A	3μm	4.6 × 100mm
HCA030U046X15022A	3μm	4.6 × 150mm
HCA030U046X25022A	3μm	4.6 × 250mm
HCA050U021X05023A	5μm	2.1 × 50mm
HCA050U021X10023A	5μm	2.1 × 100mm
HCA050U021X15023A	5μm	2.1 × 150mm
HCA050U021X25023A	5μm	2.1 × 250mm
HCA050U046X05023A	5μm	4.6 × 50mm
HCA050U046X10023A	5μm	4.6 × 100mm
HCA050U046X15023A	5μm	4.6 × 150mm
HCA050U046X25023A	5μm	4.6 × 250mm
HCA050U046X25045A	5μm	4.6 × 250mm

Note: HCA050U046X25045A is the original SCX column.