

Electrophoretic separation of proteins

Separation into polyacrylamide gel according to molecular weight. In order to separate the proteins of lower molecular weight, use of more concentrated gel is required.


Transfer of proteins

Transfer from gel onto membrane followed by:

- Blocking;
- Applying a primary antibody specific for your protein of interest;
- Applying secondary antibody that will recognize the primary antibody.

Role of protein binding

Set up for transfer



TYPICALLY TWO TYPES OF MEMBRANES:


- Nitrocellulose (hydrophilic)
- Polyvinylidene fluoride (hydrophobic)

Protein to membrane binding interactions:
hydrophobic, electrostatic, dipolar



Detection of proteins

Proteins can be detected by immunodetection methods which use enzyme conjugated/labeled secondary antibodies. When the enzyme substrate is added, a product is formed. This product can be detected by fluorescence, colorimetrically, or by chemiluminescence. Enhanced chemiluminescence (ECL) produces light as a by-product when the substrate is catalyzed by the enzyme. This light is then captured on X-ray film or by a digital imaging system.



Transfer Membranes

Unlike many transfer membrane suppliers, GVS Filter Technology manufactures all of its nitrocellulose, PVDF membranes that we sell. Our transfer membranes are used in key research and testing done around the world.

PVDF

The PVDF-Plus is a naturally hydrophobic transfer membrane designed to deliver the highest binding capacity and lowest background in protein analysis applications.

Nitrocellulose

GVS Filter Technology manufactures two pure nitrocellulose transfer membranes. NitroBind is the classic unsupported pure nitrocellulose membrane used for all protein and immunoblotting applications. NitroPure is a supported pure nitrocellulose membrane combining the characteristics of nitrocellulose with the strength of nylon. It outperforms standard nitrocellulose in reprobing applications of DNA/RNA/Protein when extensive handling is required. Most other suppliers buy membranes from a manufacturer and sell them under their label. The chart below shows which of the market leaders actually make the transfer membranes they sell.

| Membrane | Features | Benefits | Choose by Detection Systems | Choose by Procedures |
|-----------|-------------------------------|--|---|---|
| NitroPure | Pure Nitrocellulose | Pure nitrocellulose is the membrane of choice for protein and immunoblotting techniques, as well as any other procedures that require optimum resolution. Binding Capacity - 100 µg/cm ² | Radiolabeled, Chromogenic and Chemiluminescent Detection Systems | Westerns Protein & Immunoblotting Northerns Southern |
| NitroPlus | Supported Pure Nitrocellulose | Supported pure nitrocellulose is used in procedures requiring the highest sensitivities, low backgrounds and rigorous handling. The membrane can be reprobed many times. Binding Capacity 100 µg/cm ² | Radiolabeled Detection Systems Chemiluminescent and Biotinylated Detection Systems | Northern Southern Multiple Rehybridizations Colony and Plaque Lifts |
| PVDF-Plus | Hydrophobic PVDF Membrane | Hydrophobic PVDF membrane is designed for protein sequencing, western transfers and amino acid analysis. Binding Capacity - 125 µg/cm ² | Chemical compatibility allows the use of all commonly used stains | Western Transfers Protein Sequencing Amino Acid Analysis |

Selecting a Membrane

Selecting the appropriate membrane is critical to the success of a nucleic acid or protein transfer procedure.

GVS Filter Technology manufactures many types of membranes for hybridization technology, each exhibiting different performance characteristics which can directly affect the outcome of a specific technique. Below are some of the more frequently performed procedures and features of hybridization membranes.

Rehybridizations

GVS Filter Technology manufactures membranes recommended for rehybridization procedures: Magna Nylon, NitroPlus Nitrocellulose supported and NitroPure, a supported pure nitrocellulose. NitroPlus Nitrocellulose supported membranes can be most frequently reprobed. On nylon membranes, the number of reprobing steps is a function of the amount of hydrolysis to which the membrane is exposed during the protocol, and the additive effects of hot water, sodium hydroxide and an acidic environment. Sodium hydroxide solutions deteriorate the nylon matrix and are not recommended in procedures where reprobing steps are required. The polyester support web used in manufacturing NitroPure allows the membrane to be reprobed several times. Because the binding capacity of nitrocellulose is less than that of nylon (100 µg/cm² vs. 400 µg/cm²), the potential number of rehybridizations is fewer. See pages 105-107 for more details.

UV Crosslinking

For covalent binding of nucleic acids to a transfer membrane, GVS Filter Technology membranes can be UV Crosslinked by following the manufacturer's instructions. It is particularly recommended when working with short fragments, small samples, or low numbers of base pairs, because of the improved resolution this technique offers.

Protein Blotting

NitroPure nitrocellulose and PVDF-Plus membranes are recommended for use in protein blotting. Nitrocellulose membranes are able to be more thoroughly blocked, reducing the high background potential associated with protein blotting. PVDF membranes are more resistant to the harsh chemicals used in Edman degradation.

Alkaline Blotting

For more rapid transfers, an alkaline blotting procedure can be used with MagnaProbe or MagnaCharge membranes. Alkaline blotting is not recommended when reprobing is required. Please see page 125 for more details.

Staining Procedures

NitroPure, NitroPlus and PVDF-Plus membranes are recommended for procedures that require a staining step with India Ink, Coomassie Blue, Colloidal Gold, or any other commonly used stain. Nylon membranes irreversibly bind many stains.

Reducing Backgrounds

There are many sources of background problems, or low signal-to-noise ratios. Some of the most common include: contaminated probes, contaminated hybridization solutions, and incorrectly chosen stringency levels. Nonfat milk should not be used as a blocking agent as it may increase nonspecific binding. GVS Filter Technology membranes are all manufactured by strict quality control procedures, ensuring a uniform membrane with consistently low backgrounds. Please refer to pages 91-93 for more details.

Troubleshooting Common Blotting Problems

Many blotting problems can be eliminated by observing the following recommendations.

Blotchy or incomplete transfers are caused by poor contact between the gel and the membrane. Even after careful smoothing of the membrane to the gel, incomplete degassing of transfer solutions can cause air pockets to form. Evolving gas from Tris or, in the case of protein transfers, methanol, can disrupt the tight contact necessary between the membrane and the gel for successful transfers. Smeared or skewed bands are often caused by uneven contact between the gel and the membrane, or the membrane and the chromatography paper. To avoid this problem, roll a pipet down the membrane after it has been applied to the gel, and once again over the chromatography paper after it has been applied to the membrane. Do not move the membrane until the transfer is complete, as this will cause smearing.

Protocols for Protein Applications

NitroPure and NitroPlus Nitrocellulose Membranes

Gel Preparation

Western (Protein) Blotting

Gels should be stained after transfer with Coomassie Blue, Fast Green, Amido Black, or any other appropriate stain.* Soak the gel for 1 hour in a transfer buffer made of: 25 mM Tris-HCl/pH 8.0, 0.15 M glycine, 20% methanol. *GVS Filter Technology does not recommend staining before transfer. Proteins may precipitate in the membrane and not be able to transfer.

Transfer Membrane Preparation

Completely soak the membrane in deionized water, and then in transfer buffer.

Electroblotting

Assemble the membrane and gel in the electroblotting unit. Place the membrane on the anode (positive) side of the gel. Transfer according to manufacturer's instructions. Remove and wash thoroughly with transfer buffer.

Capillary Blotting

Prepare gel assembly by the method of Southern (see page 115). Transfer for 2 hours to overnight. Use transfer buffer of 10 mM Tris-HCl/pH 7.5. After the transfer step, determine transfer efficiency by staining the blot or gel by standard methods.

Blocking Procedures

Step 1: First Wash

Block the blot in PBS buffer (0.9% NaCl, 10 mM sodium phosphate/pH 7.2) containing 5% BSA, Tween 20 or high purity gelatin for 1 hour, with gentle agitation.

Step 2: Primary Antibody Binding

Remove the PBS buffer solution from blot completely. Dilute the first antibody in 50 ml of fresh PBS buffer solution. Incubate the blot in the PBS blocking buffer/antibody solution for 1 hour at 37°C with gentle agitation. Use a ratio of 5-10 ml of solution to 100 cm² of membrane.

Step 3: Second Wash

Wash the membrane in 100 ml of fresh PBS buffer solution (without antibody) with 0.1-0.3% Tween-20. Agitate in a shaker for 5 minutes. Repeat the wash step 2 times. (Note: Increasing the number of short washes reduces the potential for high backgrounds).

Detection

Thoroughly remove the PBS buffer solution and overlay the blot with an antiserum (second) antibody, or with protein A (radiolabeled or enzyme linked) for 1-2 hours at room temperature with gentle agitation. The final concentration of radiolabeled second antibody solution should be $1-2 \times 10^5$ dpm/ml of PBS buffer solution. Enzyme-linked second antibody solutions should be made at a 1:1000 titer in PBS buffer solution. Repeat the wash step described in the procedure above.

Signal Development

The choice of signal development method is dependent on the type of probe used. Radiolabeled probes are developed and quantitated by autoradiography. Enzyme-conjugated labels (horseradish peroxidase or alkaline phosphatase) are developed and quantitated with the appropriate substrate solution.

Probe Removal (Stripping)

Do not allow the filter to become dry, or irreversible binding of the probe will result.

Wash the membrane at 60°C for 30 minutes in 0.05 M sodium phosphate/pH 6.5, 10.0 M urea, 0.1 M 2-mercaptoethanol, or wash the membrane in 0.2 M glycine-HCl, 0.5 M NaCl for 5 minutes. Rinse in 0.1 M NaOH or 0.5 M Tris for 10 minutes.

PVDF-Plus

Western Blot General Protocol

Main Solutions and Reagents for running; transfer and blocking

Running buffer 10X:

- ▲ Tris base: 250 mM
- ▲ Glycine: 1.90 M
- ▲ SDS: 1%.

The pH of the buffer should be 8.3 and no pH adjustment is required. Store the running buffer at room temperature and dilute before use.

Running buffer 1X:

- ▲ 10% 10X Running buffer
- ▲ 90% DW H₂O

Tris Glycine Buffer 1X:

- ▲ 25 mM Tris base
- ▲ 190 mM Glycine

Transfer Buffer:

- ▲ 20% MeOH
- ▲ 0.25X Tris Glycine buffer

Phosphate Buffered Saline (PBS) 1X:

- ▲ 137 mM NaCl
- ▲ 2.7 mM KCl
- ▲ 10 mM Na₂HPO₄
- ▲ 1.8 mM KH₂PO₄

PBS Tween (PBST) 1X:

- ▲ 0.05% Tween
- ▲ 99.95% PBS 1X

Blocking Buffer:

- ▲ 5% skim milk (or Bovine Serum Albumin - BSA)
- ▲ 95% PBST

PROCEDURE

Electrophoresis – protein separation

1. Prepare appropriate SDS-Polyacrylamide (SDS-PAGE) gel for electrophoresis.
Type of SDS-PAGE gels according to the protein size; the lower is the protein size, the higher concentration of gel should be used.
2. Prepare the sample to be loaded in the wells of SDS-PAGE gel.
Preparation of the sample and the sample buffer depends on the type of the protein and manufacturer's recommendations.
3. Load protein marker and equal volumes of protein sample into corresponding wells of SDS-PAGE gel.
Fill the empty wells with the sample buffer.
4. Fill the electrophoresis tank with running buffer.
5. Run the gel in following conditions:
 - a. 120 V for 20-30 minutes (or until the sample reaches the stacking gel);
 - b. 180 V for 30-45 minutes (separation of the proteins under constant voltage).

Electrotransfer of Proteins

1. In case of PVDF membrane perform membrane equilibration by:
 - a. Immersing membrane in Methanol for 1 minute;
 - b. Followed by immersion of membrane in DW water for 5 minutes;
 - c. Followed by immersion of membrane in Transfer Buffer for 10 minutes.*Membrane must be wet at all times.*
2. Assemble the transfer sandwich according to scheme presented in Figure 1.
Ensure there are no bubbles between the gel and the membrane.
3. Place the cassette in the transfer tank and fill the Electroblotting tank with the transfer buffer (ensure that the sandwich is covered with the buffer).
4. Run the Electroblotting for 1 hour at 120 V in an ice bath.
Running conditions might need optimization.

Blocking and antibody incubation

1. Incubate membrane for 1h in the blocking buffer at room temperature or overnight at 4°C with constant agitation.
The active side of the membrane must always be in contact with the solution.
2. Place the blot in the primary antibody solution and incubate with agitation for 1 hour at room temperature.
The solution should move freely across the surface of the membrane (dilution of the antibody depends on the producer recommendation).
3. Wash membrane by:
 - a. Immersion in PBS-Tween (PBST) for 10 minutes with agitation;
 - b. Immersion in PBS-Tween (PBST) for 5 minutes with agitation (2 times).
4. Place the blot in the secondary antibody solution (HRP conjugates) and incubate with agitation for 45 minutes at room temperature.
Dilution of the antibody depends on the producer recommendation.
5. Wash the membrane according to the washing steps described in point 3 of Blocking and antibody incubation section.

Detection via chemiluminescence

1. Prepare a 1:1 mixture of chemiluminescent substrate (ECL HRP, depending on sensitivity choose Light Wave; Light Wave Plus or Light Wave Max) .
2. Place the blot in the container with substrate and incubate for 3 minutes.
3. Remove the excess of the solution off the membrane.
4. Place membrane in blot development folder and gently smooth out all the bubbles using a roller.
5. Expose the film to the imaging system.



Troubleshooting Guide and Application Tips

Problems and Solutions

Unsuccessful Rehybridizations

My membrane is deteriorating during the rehybridization procedure?

If so, what type of membrane are you using? GVS Filter Technology manufactures one type of membrane recommended for rehybridization procedures: NitroPlus, a supported pure nitrocellulose. A more resilient membrane during applications requiring multiple reprobings. Nitropure (a supported nitrocellulose) was developed for this reason. The polyester support web used in manufacturing NitroPlus allows the membrane to be reprobed several times.

My application demands an extensive number of reprobings and I'm losing signal?

If so, what type of membrane are you using? Because the binding capacity of nitrocellulose is less than that of nylon (100 $\mu\text{m}/\text{cm}^2$ vs. 400 $\mu\text{m}/\text{cm}^2$), the potential number of rehybridizations is fewer as compared to nylon membranes. The number of reprobing steps is a function of the amount of hydrolysis to which the membrane is exposed during the protocol, and the additive effects of hot water, sodium hydroxide and an acidic environment.

My probe is not stripping from the membrane, how should I change my procedure?

Did you let the membrane dry after the initial probe was applied? Drying causes irreversible binding of DNA to microporous membranes. If this has occurred, look through the helpful tips listed below.

My probe won't strip from the membrane, how can I rescue this blot?

Try preparing a new probe and using a different detection protocol. For example, if you prepared a biotinylated probe and detected with a streptavidin conjugate, omit the biotin-streptavidin step during rehybridization by using a directly conjugated probe, such as an alkaline phosphatase conjugated probe. If you used a radioactive probe, use a chemiluminescent system to detect after the next hybridization (or vice versa). If you have enough time and are using radioactive probes (e.g., pgs 106-107), simply let your first probe decay before the second round of hybridization.

Signal Problems

The nucleic acid did not transfer completely to the membrane, what should I do?

Blotchy or incomplete transfers are caused by poor contact between the gel and the membrane. Even after careful smoothing of the membrane to the gel, incomplete degassing of transfer solutions can cause air pockets to form. Evolving gas from Tris or, in the case of protein transfers, methanol, can disrupt the tight contact necessary between the membrane and the gel for successful transfers.

My signal is low, what are the common reasons for this?

When you have low signal, it is best to check your reagents by performing extra controls. The most common reason for poor signal is a bad probe. Prepare a new probe and perform a dot blot comparing the old and new probes. Do you see a difference between the probes? Even nonradioactive probes can deteriorate during storage. Is the signal weak for the new probe as well? Then your detection enzymes may be bad or the reagents used to prepare the probe are bad. You might also blot a small amount of unlabeled complementary DNA and hybridize to the new probe. Are you seeing signal from the blotted probe but not the hybridized DNA? If so there could be a problem with your hybridization protocol, such as the wash temperature or your buffers. If you're using nonradioactive detection methods, test your enzymes and substrates as well.

Background Problems

Everything was working fine and now suddenly I have high backgrounds, Why?

Did you make up a new probe? If so, was there adequate separation of the unincorporated label from the incorporated? Are you using old solutions? There may be contamination. Usually in these cases it is best to prepare new solutions, new probes and use new reagents. This is often the fastest way to get your system working again.

Miscellaneous

My membrane changed color during my blotting procedure, should I be concerned?

No. Slight color changes in GVS Filter Technology new positively charged membranes are expected and have no effect on results. These color changes will vary according to the blotting procedure used and the pH of solutions. GVS Filter Technology uses this color change to ensure quality during the manufacturing procedure.

High Backgrounds

Poor agitation during prehybridization and hybridization steps can lead to insufficient blocking of the entire membrane. Due to the strength of the internal support web, NitroPlus can withstand higher levels of agitation without tearing or ripping. Incorrect probe concentration can occur when using dextran sulfate in hybridization or prehybridization solutions. Dextran sulfate causes the effective concentration of the probe to increase because it excludes the probe from the volume of solution the dextran sulfate polymer occupies. When using dextran sulfate, lower the probe solutions to less than 10 ng/ml of the solution. When not using dextran sulfate, maintain the optimum probe concentration at 25-40 ng/ml of solution. Residual agarose on membranes can cause a fuzzy background to appear on blots. Be sure to wash nylon membranes with 5 x SSPE at 60°C, after the immobilization step. Due to the strength of the membrane, supported membranes (NitroPlus) can be more easily washed without tearing or ripping.

Troubleshooting Gel Casting Procedures

Troubleshooting blotting problems begins with the correct gel casting procedures. Skewed, streaked, incomplete, or nonuniform transfers can be the results of poorly cast gels. The following recommendations are made for setting up the gel.

Gels greater than 4mm thick can interfere with the free transfer of nucleic acids.

Be sure that the gel tray is level before casting the gel. If the surface is not level, non-uniform transfers may result. Maintain a gel casting temperature of 55-70°C degrees, and be sure that the gel particles are completely dissolved. Undissolved agarose particles can result in streaked or skewed bands. Immediately after gel casting and solidification, submerge the gel slab in electrophoresis buffer. This will prevent the formation of an impermeable "skin" over the surface of the gel which can inhibit transfer of nucleic acids from the gel. After setting up the blotting assembly, be sure to:

- Invert the gel so that the underside of the gel is the side in contact with the membrane.
- Allow the transfer solutions enough time to "breathe," so that they may degas completely. Incompletely degassed transfersolutions evolve gas after the blotting assembly is set up, and can cause air bubbles between the membrane and gel that can impede the transfer of nucleic acids.

Probe Related Background Problems

While there are several ways to decontaminate probe solutions, the following methods are two of the most efficient. The second method can be rapidly performed with minimum effort.

Method 1: Phenol/Chloroform extract the probe to remove unincorporated nucleotides, proteins, and other contaminants.

Method 2: Clean the probe by adding a small volume of the hybridization buffer to the probe and filtering it through an ABluo 25AS low protein binding cellulose acetate syringe filter. Contaminants in the probe solution will be held back by the 0.2µm filter with no probe loss caused by nonspecific binding to the filtration membrane. Probe length is also a factor contributing to background levels seen on transfer membranes. Between 250-800 base pairs is the recommended optimum length of a probe; probe lengths smaller or larger than this can lead to a low signal-to-noise ratio. Probes smaller than 250 base pairs often bind poorly and may require less stringent hybridization and wash procedures. Probes larger than 800 base pairs may contain a wider variety of size classes, which can lead to extraneous binding to the transfer membrane.

Hybridization Solution Related Background Problems

Contaminated hybridization solutions are another common source of background problems. Hybridization solutions should be filtered with a pure cellulose acetate ABluo 25AS syringe filter, to remove contaminants.

Additionally, all solutions and buffers should be made fresh before each transfer with sterile, double-distilled, deionized water, and very high grade reagents. After fresh buffers are made, they should be filtered with an ABluo 25AS syringe filter to ensure that no contaminants remain in the solution. Formamide-based hybridization solutions are a frequent source of background noise, and the formamide must be freshly made and deionized.


Optimized Blocking Solutions

A concentration of 5-7 x Denhardt's solution is recommended for use with nylon membranes. Exceeding this level can lead to quenching of the signal.

Backgrounds Associated with Reprobing

A follow-up autoradiograph after probe removal is strongly recommended to determine if the probe has been fully stripped. Otherwise, backgrounds can appear in blots that have not been fully erased.

Nitrocellulose



GVS Nitrocellulose Pure Transfer Membrane is the membrane of choice for all protein or immunoblotting applications. The high sensitivity of GVS Nitrocellulose Transfer Membrane ensures excellent results in all transfers, especially in protein blotting. Unlike PVDF, nitrocellulose wets out naturally, does not require methanol, and will not turn hydrophobic during the transfer process.

Nitrocellulose is very easily blocked and does not need the many blocking steps required with PVDF.

Excellent results will be obtained with all detection systems: antibody/antigen, radiolabeled, biotinylated, and chemiluminescent, giving you a great amount of flexibility in designing your procedure.


Supplied in various porosity and format.

Features & Benefits

- ▲ For procedures that require optimum resolution
- ▲ Membrane of choice for protein or immunoblotting applications
- ▲ Low background, easily blocked
- ▲ BSA binding capacity up to 100 µg/cm²
- ▲ Wets out naturally
- ▲ Compatible with all detection systems

Typical Applications

- ▲ Western Blotting
- ▲ Protein & immunoblotting
- ▲ Northern Blotting
- ▲ Southern Blotting
- ▲ Dot/slot blotting
- ▲ Radiographic, chromogenic and chemiluminescent detection systems

| Product | Competitors |
|---|--|
|  NITRO Pure Unsupported Nitrocellulose | Amersham HyBond-C - BioRad Nitrocellulose - Millipore Immobilon-NC Plus - Schleicher & Shuell (S&S) Protran |

Ordering information

| | Dimensions (mm) Packaging | 70x84 mm 10/pk | 100x100 mm 10/pk | 150x150 mm 5/pk | 200x200 mm 25/pk | 200x3000 mm 1/pk | 300x3000 mm 1/pk |
|------------|------------------------------|-------------------|---------------------|--------------------|---------------------|---------------------|---------------------|
| Pore sizes | 0.22 µm | 1213991 | 1213999 | 1215463 | 1215392 | 1215469 | 1215458 |
| | 0.45 µm | 1213888 | 1213314 | 1215476 | 1221976 | 1215483 | 1215471 |

Supported Nitrocellulose




GVS Supported Nitrocellulose Transfer Membrane combines the binding characteristics of nitrocellulose membrane with the strength of nylon membrane. It can be easily used in any protocol utilizing unsupported nitrocellulose transfer membrane. Supplied in various porosity and format

Features & Benefits

- Supported for procedures requiring rigorous handling
- Strong - will not curl, bend or crack after baking
- High sensitivities, low backgrounds
- Multiple reprobing
- BSA binding capacity up to 100 µg/cm²
- Triton Free

Typical Applications

- Northern Blotting
- Southern Blotting
- Multiple re-hybridizations
- Colony/plaque lifts
- Dot/slot blotting
- Radiographic detection systems
- Chemiluminescent detection systems
- Biotinylated detection systems



All lanes : Anti-Furin antibody [EPR14674] (ab183495) at 1/5000 dilution

Lane 1 : HepG2 whole cell lysate

Lane 2 : HeLa whole cell lysate

Lane 3 : U87-MG whole cell lysate

Lane 4 : Caco-2 whole cell lysate

Lysates/proteins at 20 µg per lane.

Secondary

Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/1000 dilution

Predicted band size : 87 kDa

| Product | Competitors |
|---------|-------------|
|---------|-------------|



Amersham HyBond-C Extra - Amersham HyBond-C Super
Biorad Supported Nitrocellulose


Ordering information

| | Dimensions (mm) Packaging | 70x84 mm 10/pk | 100x100 mm 10/pk | 150x150 mm 5/pk | 200x200 mm 5 /pk | 200x3000 mm 1/pk | 300x3000 mm 1/pk |
|------------|------------------------------|-------------------|---------------------|--------------------|---------------------|---------------------|---------------------|
| Pore sizes | 0.22 µm | | 1214560 | 1212669 | 1212689 | 1212690 | 1212632 |
| | 0.45 µm | 1214978 | 1213943 | 1212596 | 1212597 | 1212602 | 1212590 |


Polyvinylidene Fluoride PVDF



GVS



Competitor



Images were obtained by following GVS Western Blot General Protocol

Cell Lane: HeLa Whole Cell

Detection substrate: Light Wave Plus

Primary antibody: Beta Actin Polyclonal Antibody (dilution 1:1000)

Secondary antibody: Goat Anti-Rabbit IgG Antibody (H+L) (dilution 1:10000)


Analyzed protein: Beta actin, MW: 42 kDa

GVS PVDF is a naturally hydrophobic, unsupported transfer membrane. It has a high binding capacity, which prevents protein from passing through the membrane, and a low background that provides for an excellent signal-noise ratio. It also has exceptional tensile strength, preventing it from cracking, tearing, breaking or curling. This membrane also has broad chemical compatibility, which is important when used with common stains such as Amido Black, Colloidal Gold, Coomassie Blue, India Ink and Ponceau-S. GVS PVDF will not degrade, distort or shrink when a high concentration of methanol is used for destaining.

Its exceptional strength, high binding capacity and chemical compatibility make GVS PVDF ideal for use in Western blotting, immunoblotting, and solid phase assays and plaque lifts.

Features & Benefits

- ▲ Superior strength: Can withstand aggressive handling or be used with automated equipment without breaking or tearing
- ▲ Low extractables: Ensures tests will be clean with consistent results
- ▲ Exceptional sensitivity: Detects low-level components
- ▲ Hydrophobic: For high protein binding
- ▲ Lot-to-lot consistency: Quality checks ensure consistent binding for dependable results every time
- ▲ BSA protein binding capacity : 125 µg/cm²
- ▲ High range of chemical: Resistant to most commonly used chemicals compatible with chemically aggressive solvents

| Product | Competitors |
|--|--|
|  PVDF Pure Unsupported PVDF | Millipore Immobilon-P - Amersham HyBond - BioRad PVDF Schleicher & Shuell (S&S) Westran |

Typical Applications

- ▲ Western blotting
- ▲ Immunoblotting
- ▲ Solid phase assays
- ▲ Amino acid or protein analyses

Ordering information

| | Dimensions (mm) | 70x84 mm | 100x100 mm | 150x150 mm | 200x200 mm | 200x3000 mm | 300x3000 mm |
|------------|-----------------|----------|------------|------------|------------|-------------|-------------|
| | Packaging | 10/pk | 10/pk | 5/pk | 5/pk | 1/pk | 1/pk |
| Pore sizes | 0.22 µm | 1214588 | | 1215037 | 1215032 | 1214726 | 1214429 |
| | 0.45 µm | 1213992 | 1212644 | 1212636 | 1212637 | 1212783 | 1212639 |